



## 6 Core PFA Insulated, Tin Plated Copper Screen, PRT Sensor Cable / Wire



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**6 Core PFA Insulated, Tin Plated Copper Screen, PRT Sensor Cable / Wire**

High temperature to 260°C, extruded PFA Teflon® construction. Commonly used for Pt100? thermometers, also suitable with other types such as Pt130, Pt500 & Pt1000?

**Specifications****Specifications**

<b>General Description</b>	PFA insulated 7/0.2mm tin plated copper conductors, twisted cores with nylon cord, tin plated copper braided screen with PFA insulated outer jacket
<b>Cable Length</b>	10m, 25m, 50m, 100m
<b>Cable Type</b>	PFA Teflon® construction
<b>Core / Strands</b>	7/0.2mm
<b>Insulation Material</b>	PFA Insulated
<b>Max. Temperature</b>	+260°C
<b>Min. Temperature</b>	-75°C
<b>Number of Wires</b>	6
<b>Standards Met</b>	Cores colour coded Red & White in accordance with IEC-751, black outer jacket