

KNE Stainless Steel Terminal Head



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Many alternative types of terminal heads are available to meet the requirements of various applications.

Variations exist in size, material, accommodation, resistance to media, resistance to fire or even explosion and in other parameters.

KNE Stainless Steel Terminal Heads provide robust protection for electrical connections in harsh environments. Their corrosion resistance, durability, and secure sealing make them suitable for a wide range of industrial, chemical, food and beverage, and marine applications.

Proper installation and regular maintenance ensure reliable performance and long-term protection of sensitive electrical components.

Specifications

Specifications

Product Code	XE-3823-001
Head Type	KNE Stainless steel Terminal Head
Material	316 Stainless Steel
Cable Entry	M20 x 1.5mm
Fixing (Centres)	33mm & 40mm
Max. Temperature	170°C
IP Rating	IP67 rated
Process Entry	1/2 " BSPP
Cap Retaining Type	Screw cap and retaining screw
Series	Terminal Heads