

## Glassfibre insulated Flat Pair IEC Exposed Junction Thermocouple - Type J with Bare Tails Termination - 1/0.315mm x 2M

## Image Coming Soon

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.



## **Specifications**

## **Specifications**

Product Code	XE-3571-001	
General Description	Good temperature resistance (not suitable in fluids)	
Sensor Type	Exposed Junction	
Thermocouple Type	J	
Cable Length	2M	
Cable Type	Glassfibre	
Core / Strands	1/0.315mm	
Termination Type	Bare Tails	
Max. Temperature	+350°C (short periods up to 400°C)	
Min. Temperature	-60°C	
Tolerance	Class 2	
Standards Met	IEC-584-3	