



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

PFA Insulated IEC Fine Gauge Exposed Junction Thermocouple (0.076mm conductors) - Type K - 0.5M Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Thermocouple type K. Fast response, welded exposed junction available with bare tails. Ideal for test & development applications.

Specifications**Specifications**

Product Code	XE-3509-001
General Description	Fast response, welded tip exposed junction
Sensor Type	Exposed Junction
Thermocouple Type	K IEC
Cable Length	0.5M
Cable Type	Teflon® PFA insulated wires
Core / Strands	0.003" (0.076mm) conductors
Termination Type	Bare Tails
Max. Temperature	260°C
Min. Temperature	-75°C
Standards Met	IEC-584-3
Series	Exposed Junctions