

## Thermo Attachment Pads in PTFE glass cloth



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Self-adhesive pads in PTFE glass cloth for locating welded tip thermocouples

Rated -73°C to 260°C

Pad size 13x18mm

Sheet contains 20 pads

## Specifications

<b>Product Code</b>	XE-3599-001
<b>Series</b>	ATTACHMENT PADS