



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Industrial Style Mineral Insulated Thermocouple with KNE Terminal Head - Type K - 6.0mm Diameter - 600mm Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Need Your Temperature Sensor calibrated?

Traceable to ISO17025 (UKAS) calibration standards, Labfacility offer a choice of a 3 or 5 point traceable calibration using our in-house calibration facility, you can select any temperature points between -10°C & +200°C.

Please be aware that if you add a calibration option to your order you will be contacted regarding delivery information.

Specifications**Specifications**

Product Code	XE-3255-001
General Description	Highly flexible, sheath can be bent/formed to suit many applications and processes
Sensor Type	Mineral insulated Thermocouple fitted with KNE terminal head
Thermocouple Type	K
Thermocouple Junction	Insulated hot junction
Probe Diameter	6.0mm
Probe Length	600mm
Sheath Material	310 stainless steel sheath
Head Type	KNE, Cast aluminium finish, screw type lid with chain
Gland	M16 x 1.5mm Plated brass gland fitted
Max. Temperature	1100°C
Min. Temperature	-40°C
IP Rating	IP67