

Industrial Style Mineral Insulated Thermocouple with KNE Terminal Head - Type J - 6.0mm Diameter - 150mm Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.



Datasheet

Need Your Temperature Sensor calibrated?

Traceable to ISO17025 (UKAS) calibration standards, Labfacility offer a choice of a 3 or 5 point traceable calibration using our in-house calibration facility, you can select any temperature points between -10°C & +200°C.

Please be aware that if you add a calibration option to your order you will be contacted regarding delivery information.

Specifications

Specifications

XE-3257-001
Highly flexible, sheath can be bent/formed to suit many applications and processes
Mineral insulated Thermocouple fitted with KNE terminal head
J
Insulated hot junction
6.0mm
150mm
310 stainless steel sheath
KNE, Cast aluminium finish, screw type lid with chain
M16 x 1.5mm Plated brass gland fitted
1100°C-40°C
IP67