


ANSI Glassfibre insulated Exposed Junction Thermocouple with Bare Tails - Type K - 1M Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.



Specifications

Specifications

Product Code	XE-3545-001
General Description	Fast response, welded exposed junction
Sensor Type	Good temperature resistance (not suitable in fluids)
Thermocouple Type	K
Thermocouple Junction	Exposed
Cable Length	1M
Cable Type	glassfibre insulated flat pair cable
Core / Strands	1/0.3mm
Cable Termination	Bare tails
Max. Temperature	+350°C (short periods up to 400°C)
Min. Temperature	-60°C
Tolerance	Class 2
Standards Met	ANSI MC96.1
otalidal do mot	74401141000.1