



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

ANSI Mineral Insulated Thermocouple - Type K with Miniature Plug - 3.0mm x 1000mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Need Your Temperature Sensor calibrated?

Traceable to ISO17025 (UKAS) calibration standards, Labfacility offer a choice of a 3 or 5 point traceable calibration using our in-house calibration facility.

Please be aware that if you add a calibration option to your order you will be contacted regarding delivery information.

Specifications**Specifications**

Product Code	XE-3018-001
Thermocouple Type	K
Thermocouple Junction	Hot junction insulated from sheath
Probe Diameter	3.0mm
Probe Length	1000mm
Sheath Material	Stainless steel 310 sheath
Cable Termination	Miniature Plug
Max. Temperature	1100°C
Standards Met	ANSI