


## General Purpose Pt100 Silicone Rubber Lead - Type PRT - 6.0mm Diameter - 50mm Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Specifications****Specifications**

<b>Product Code</b>	XE-3630-001
<b>General Description</b>	Suitable for general purpose industrial & commercial use
<b>Sensor Type</b>	Pt100 (100 Ohms @ 0°C) to IEC 751, Class B, 4 wire
<b>Probe Diameter</b>	6.0mm
<b>Probe Length</b>	50mm
<b>Sheath Material</b>	316 stainless steel sheath
<b>Cable Length</b>	2M
<b>Cable Type</b>	Flexible silicone rubber insulated
<b>Core / Strands</b>	7/0.2mm
<b>Cable Termination</b>	Bare tails
<b>Termination Type</b>	4-wire tails allows connection to any PT100 instrument 2, 3 or 4 wire
<b>Max. Temperature</b>	200°C
<b>Min. Temperature</b>	-50°C
<b>Accuracy</b>	Class B
<b>Number of Wires</b>	4
<b>Series</b>	General Purpose Pt100