



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Stainless Steel Sheathed Pt100 Sensor - Type PRT - 1M Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-3665-001
General Description	Sheathed Pt100 element in a thin, flat stainless steel body
Sensor Type	Quick thermal response while element remains protected
Element	Pt100
Cable Length	1M
Cable Type	Teflon® insulated twin twisted lead
Core / Strands	7/0.2mm
Cable Termination	Bare tails
Max. Temperature	250°C
Min. Temperature	-50°C
Accuracy	Class B
Number of Wires	2
Series	Sheathed Pt100