


Type K IEC Right Angled - Copper Disc Tip



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Stainless steel sheath with a moulded handle and 2 metres of coiled extension cable terminated in a miniature thermocouple plug.

Specifications

Product Code	XE-3817-001
General Description	Surface temperature spring loaded thermocouple with right angled probe and copper disc tip
Thermocouple Type	K
Calibration	IEC
Cable Length	2m
Cable Termination	Fitted Plug
Max. Temperature	+600°C
Dimensions	Probe Ø8mm(17mmL) Tip x 60mm Hot Leg and Ø6mm x 180mm Cold Leg with Ø4.5mm copper disc.
Series	Handheld Probe