


## Type K IEC Right Angled Surface



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Stainless steel sheath with a moulded handle and 2 metres of coiled extension cable terminated in a miniature thermocouple plug.

## Specifications

<b>Product Code</b>	XE-3819-001
<b>General Description</b>	Surface temperature fast response thermocouple
<b>Sensor Type</b>	Right angled Probe
<b>Thermocouple Type</b>	K
<b>Calibration</b>	IEC
<b>Cable Length</b>	2m
<b>Cable Termination</b>	Fitted Plug
<b>Max. Temperature</b>	+900°C
<b>Dimensions</b>	Probe Ø12mm Tip x 40mm Hot Leg and Ø6mm x 180mm Cold Leg. Coiled element diameter 7mm
<b>Tip</b>	Ceramic Tip, Coiled Element
<b>Series</b>	Handheld Probe