


Type K IEC Right Angled Surface



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Stainless steel sheath with a moulded handle and 2 metres of coiled extension cable terminated in a miniature thermocouple plug.

Specifications

Product Code	XE-3819-001
General Description	Surface temperature fast response thermocouple
Sensor Type	Right angled Probe
Thermocouple Type	K
Calibration	IEC
Cable Length	2m
Cable Termination	Fitted Plug
Max. Temperature	+900°C
Dimensions	Probe Ø12mm Tip x 40mm Hot Leg and Ø6mm x 180mm Cold Leg. Coiled element diameter 7mm
Tip	Ceramic Tip, Coiled Element
Series	Handheld Probe