




## PCB Flat Film Detectors (100 Ohm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

<b>Product Code</b>	XE-4480-001
<b>General Description</b>	PCB Mounting with solder pads
<b>Sensor Type</b>	Pt100 to IEC751
<b>Max. Temperature</b>	+150°C
<b>Min. Temperature</b>	-50°C
<b>Detector Dimensions</b>	3.2mm x 1.6mm, 0.9mm thickness
<b>Features #1</b>	Resists vibration & shock
<b>Features #2</b>	Thin film construction
<b>Series</b>	RTD Detectors