




PCB Flat Film Detectors (100 Ohm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

Product Code	XE-4480-001
General Description	PCB Mounting with solder pads
Sensor Type	Pt100 to IEC751
Max. Temperature	+150°C
Min. Temperature	-50°C
Detector Dimensions	3.2mm x 1.6mm, 0.9mm thickness
Features #1	Resists vibration & shock
Features #2	Thin film construction
Series	RTD Detectors