

## 100 Ohm Wire Wound Detectors - Class A - 2.8mm x 15mm



100 OHM WIRE WOUND DETECTOR



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Please note: Class 1/10 DIN lead times will vary - please enquire or you will be contacted with a lead time after placing your order.

### Specifications

## Specifications

|                            |                                                              |
|----------------------------|--------------------------------------------------------------|
| <b>Product Code</b>        | XE-4746-001                                                  |
| <b>General Description</b> | Wirewound construction for optimum performance and stability |
| <b>Sensor Type</b>         | Pt100 Single Element                                         |
| <b>Max. Temperature</b>    | 650°C                                                        |
| <b>Min. Temperature</b>    | -200°C                                                       |
| <b>Tolerance</b>           | Class A                                                      |
| <b>Dimensions</b>          | 2.8mm Width x 15mm Length                                    |
| <b>Series</b>              | RTD Detectors                                                |