

100 Ohm Wire Wound Detectors - Class A - 2.8mm x 25mm



100 OHM WIRE WOUND DETECTOR



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Please note: Class 1/10 DIN lead times will vary - please enquire or you will be contacted with a lead time after placing your order.

Specifications

Specifications

Product Code	XE-4747-001
General Description	Wirewound construction for optimum performance and stability
Sensor Type	Pt100 Single Element
Max. Temperature	650°C
Min. Temperature	-200°C
Tolerance	Class A
Dimensions	2.8mm Width x 25mm Length
Series	RTD Detectors