




LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

100 Ohm Wire Wound Detectors - Class 1/10 DIN - 1.5mm x 8mm



100 OHM WIRE WOUND DETECTOR



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Please note: Class 1/10 DIN lead times will vary - please enquire so we can advise before placing your order.

Specifications

Specifications

Product Code	XE-4788-001
General Description	Wirewound construction for optimum performance and stability
Sensor Type	Pt100 Single Element
Max. Temperature	650°C
Min. Temperature	-200°C
Tolerance	1/10DIN relative to IEC 60751 over the range 0°C to 150°C
Dimensions	1.5mm Width x 8mm Length