

100 Ohm Dual Element Wire Wound Detectors - Class B Duplex - 2.8mm x 25mm





Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.



Specifications

Specifications

XE-4817-001
Wirewound construction for optimum performance and stability
Dual element Pt100
650°C
-200°C
Class B to IEC 60751
2.8mm Width x 25mm Length
RTD Detectors