




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

100 Ohm Dual Element Wire Wound Detectors - Class A Duplex - 2.8mm x 15mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-4846-001
General Description	Wirewound construction for optimum performance and stability
Sensor Type	Dual element Pt100
Max. Temperature	650°C
Min. Temperature	-200°C
Tolerance	Class A to IEC 60751
Dimensions	2.8mm Width x 15mm Length