

Pt100 4 wire class B Resistance Thermometer, KNE Head



[Click here to view product](#)

| | |
|-------------------------|---|
| Product Code | XE-3652 |
| Sensor Type | Pt100 (100 Ohms @ 0°C) to IEC 751, Class B, 4 wire |
| Probe Diameter | 6mm |
| Probe Length | 200mm, 400mm or 600mm |
| Sheath Material | 316 stainless steel |
| Head Type | KNE |
| Termination Type | 4 wire connection block, M20 x 1.5mm cable entry (gland included) |
| Max. Temperature | +450°C (connection head @ 170°C) |
| Min. Temperature | -50°C |
| IP Rating | IP67 |

| | |
|------------------------|---------|
| Accuracy | Class B |
| Number of Wires | 4 |

