


General Purpose Thermocouple Probe, Glassfibre stainless steel overbraid - Type K,J



General Description	Rugged thermocouple with stainless steel sheath
Thermocouple Type	K or J
Thermocouple Junction	Grounded hot junction for fast response
Probe Diameter	4.5mm
Probe Length	150mm
Cable Length	2m
Cable Type	Glassfibre insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Termination Type	Bare tails or Standard Plug
Max. Temperature	+350°C
Min. Temperature	-60°C