



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

IEC Adjustable Bayonet Thermocouple, Glassfibre Stainless Steel Overbraided Cable - Type J - Bare Tails Termination - 2M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-3405-001
General Description	Thermocouple with adjustable insertion depth
Sensor Type	6.0mm diameter tip
Thermocouple Type	J
Thermocouple Junction	Grounded tip for fast response
Cable Length	2M
Cable Type	Glassfibre insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Termination Type	Bare Tails
Max. Temperature	350°C
Min. Temperature	-60°C
Standards Met	IEC-584-3
Bayonet cap	11.5mm I.D. single slot
Thread	M6 x 1mm thread
Spring	200mm