

**IEC Pipe Surface Jubilee Clip Style  
Thermocouple, Glassfibre stainless steel  
overbraid - Type K - to fit 13-20mm Pipe Diameter**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Specifications****Specifications**

<b>Product Code</b>	XE-3914-001
<b>General Description</b>	Stainless steel pipe clamp sensor, jubilee clip style, with adjustable diameter
<b>Thermocouple Type</b>	K
<b>Thermocouple Junction</b>	Grounded tip for fast response
<b>Cable Length</b>	2M
<b>Cable Type</b>	Glassfibre stainless steel braided extension cable
<b>Core / Strands</b>	7/0.2mm
<b>Termination Type</b>	Bare Tails
<b>Max. Temperature</b>	350°C
<b>Standards Met</b>	IEC
<b>To Fit Pipe Diameters</b>	13-20mm