


## IEC Autoclave Thermocouple - Duplex Type T



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

This Steriliser Load Probe offers a Duplex Type T sensor, and is suitable for use in harsh conditions, including most steam sterilisers.

## Specifications

<b>Product Code</b>	XE-9510-001
<b>General Description</b>	First 2.5m section in a flexible stainless steel conduit
<b>Thermocouple Type</b>	Type T Duplex
<b>Thermocouple Junction</b>	2 x insulated tip junctions
<b>Probe Diameter</b>	4mm
<b>Probe Length</b>	250mm
<b>Sheath Material</b>	Stainless steel (316) sheath
<b>Cable Type</b>	Cable is silicon rubber throughout
<b>Max. Temperature</b>	+200°C
<b>Min. Temperature</b>	-50°C