




**LABFACILITY**  
TEMPERATURE & PROCESS TECHNOLOGY

## **IEC Air Probe Style Thermocouple Probe with Moulded on Miniature Plug - Type K - 110mm**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

316 Stainless Steel Air Probe with a vented sheath for general air temperature measurement.

### Specifications

## Specifications

<b>Product Code</b>	XE-0013-001
<b>General Description</b>	Vented sheath for general air temperature measurement
<b>Sensor Type</b>	Mineral insulated
<b>Thermocouple Type</b>	K
<b>Probe Diameter</b>	4.0mm
<b>Probe Length</b>	110mm
<b>Termination Type</b>	Miniature moulded flat pin plug, colour coded in accordance IEC specification
<b>Max. Temperature</b>	+250°C (Moulded-on plug +110°C)
<b>Standards Met</b>	IEC
<b>Connector Type</b>	Plug
<b>Connector Size</b>	Miniature (moulded on)
<b>Probe Material</b>	vented sheath, 316 stainless steel