

Black ABS Equipment Case



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Black ABS case comes with complete foam set which can be plucked out in 10mm cubes and modified to the users requirements.

Dimensions

Internal (LxWxD)
420 x 298 x 131mm

External (LxWxD)
432 x 363 x 138mm

Weight (Kg)
1.4

Specifications

| | |
|---------------------|-------------|
| Product Code | XE-6050-001 |
|---------------------|-------------|