

## ANSI 12 Way Thermocouple Selector Switch, Bench Mounting



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

A 12 way selector in a moulded housing with aluminium front and rear panels. The rear inputs are switched via a front mounted rotary selector switch to the single output socket.

Internally wired with thermocouple compensating cable and available with thermocouple types K, J, T, N & R/S.

**SPECIFICATION:**

Construction - Portable / Bench mounting unit in moulded housing with aluminium front and rear panels

Size - 130mm wide x 60mm high x 145mm deep

Input connectors – 12 miniature compensated thermocouple sockets

Output connectors - Single miniature compensated thermocouple socket

Available in ANSI type K, J, T, N & R/S

## **Specifications**