



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Certificate of Conformity

ISO9001

CERTIFICATE OF CONFORMITY



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

A document certified by a competent authority that the supplied good or service meets the required specifications.

Also called certificate of conformance or certificate of compliance.

Specifications

Product Code	PO-COC
--------------	--------