





LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

Ceramic Sheathed Pt100 Inserts - 4-Wire with 450mm Lead Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-3667-001
General Description	Pt100 thin film element housed in ceramic sheath
Element	Pt100
Cable Length	450mm
Cable Type	PTFE insulated wires
Core / Strands	7/0.15mm
Cable Termination	Bare Tails
Max. Temperature	250°C
Min. Temperature	-50°C
Tolerance	Class B
Detector Dimensions	Ceramic 5mm dia x 35mm long
Number of Wires	4
Series	Pt100 Inserts