






LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

IEC Fabricated Grounded Thermocouple in Stainless Steel Tube, Glassfibre stainless steel Overbraided cable - Type K,J



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

Sensor Type	316 stainless steel tube in various sizes
Thermocouple Type	K or J
Thermocouple Junction	grounded at tip for fast response
Cable Length	2m
Cable Type	Glassfibre stainless steel overbraided cable
Cable Termination	Bare wire tails
Max. Temperature	+350°C
Min. Temperature	-60°C