





# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## **IEC Fabricated Grounded Thermocouple in Stainless Steel Tube, Glassfibre stainless steel Overbraided cable - Type K - 4.0mm Diameter - 13mm Length**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Specifications****Specifications**

|                              |  |
|------------------------------|--|
| <b>Product Code</b>          | XE-5544-001                                  |
| <b>General Description</b>   | 316 stainless steel tube in various sizes    |
| <b>Thermocouple Type</b>     | K  |
| <b>Thermocouple Junction</b> | grounded at tip for fast response            |
| <b>Cable Length</b>          | 2M   |
| <b>Cable Type</b>            | Glassfibre stainless steel overbraided cable |
| <b>Cable Termination</b>     | Bare wire tails                              |
| <b>Max. Temperature</b>      | 350°C  |
| <b>Min. Temperature</b>      | -60°C  |
| <b>Standards Met</b>         | IEC  |