

**IEC Fabricated Grounded Thermocouple in  
Stainless Steel Tube, Glassfibre stainless steel  
Overbraided cable - Type J - 4.76mm Diameter -  
75mm Length**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Specifications****Specifications**

<b>Product Code</b>	XE-5550-001
<b>General Description</b>	316 stainless steel tube in various sizes
<b>Thermocouple Type</b>	J
<b>Thermocouple Junction</b>	grounded at tip for fast response
<b>Cable Length</b>	2M
<b>Cable Type</b>	Glassfibre stainless steel overbraided cable
<b>Cable Termination</b>	Bare wire tails
<b>Max. Temperature</b>	350°C
<b>Min. Temperature</b>	-60°C
<b>Standards Met</b>	IEC