


# Hygienic Pt100 Resistance Thermometer, 1.5" RJT Style Fitting



<b>General Description</b>	Ideal for hygienic processes including dairy, brewery and pharmaceutical
<b>Sensor Type</b>	Pt100, 3 or 4 wire single element class A
<b>Head Type</b>	KNEF
<b>Connection Type</b>	Pt100, 3 or 4 wire single element class A
<b>Transmitter</b>	Optional SEM-203p transmitter, ranged -50°C to +150°C
<b>Support Dia (D2)</b>	8.0mm
<b>Support Length (L2)</b>	50mm