



LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

Pt100 4 wire class B Resistance Thermometer, DIN B Head with 200mm Probe Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-5516-001
Sensor Type	Pt100 (100 Ohms @ 0°C) to IEC 751, Class B, 4 wire
Process Connection	Stainless steel plug (no thread)
Probe Diameter	6.0mm
Probe Length	200mm
Sheath Material	316 Stainless Steel
Head Type	DIN B
Max. Temperature	+450°C (connection head @ 170°C)
Min. Temperature	-50°C
IP Rating	IP67
Accuracy	Class B
Number of Wires	4
Series	Pt100 DIN B