

## 4-20mA Remote Wall Mounted Housing, Pt100 input 1 metre lead RTD / PRT



[Click here to view product](#)

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | XE-5520   |
| <b>General Description</b> | Compatible with Pt100 platinum resistance thermometers                |
| <b>Sensor Type</b>         | Surface Mounting for remote installtion                               |
| <b>Cable Length</b>        | 1m  |
| <b>Core / Strands</b>      | 7/0.2mm   |
| <b>Cable Termination</b>   | Bare tails  |
| <b>Thread</b>              | Entry Threads M16   |
| <b>Transmitter</b>         | Yes   |
| <b>Range</b>               | -50°C to +150°C or 0°C to +400°C (either can be re-ranged if desired) |
| <b>Output</b>              | 4-20 mA output  |