


## **Glassfibre insulated Flat Pair IEC Exposed Junction Thermocouple - Type K with Miniature Plug Termination - 1/0.508mm x 2M**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Need Your Temperature Sensor calibrated?**

Traceable to ISO17025 (UKAS) calibration standards, Labfacility offer a choice of a 3 or 5 point traceable calibration using our in-house calibration facility, you can select any temperature points between -10°C & +200°C.

**Please be aware that if you add a calibration option to your order, you will be contacted regarding delivery information**

**Specifications****Specifications**

<b>Product Code</b>	XE-5569-001
<b>General Description</b>	Good temperature resistance (not suitable in fluids)
<b>Sensor Type</b>	Exposed Junction
<b>Thermocouple Type</b>	K
<b>Cable Length</b>	2M
<b>Cable Type</b>	Glassfibre
<b>Core / Strands</b>	1/0.508mm
<b>Termination Type</b>	Miniature flat pin plug, colour coded in accordance IEC 584
<b>Max. Temperature</b>	+350°C (short periods up to 400°C)
<b>Min. Temperature</b>	-60°C
<b>Tolerance</b>	Class 2
<b>Standards Met</b>	IEC-584-3
<b>Connector Type</b>	Plug
<b>Connector Size</b>	Miniature
<b>Colour</b>	Green



Series	Exposed Junctions
--------	-------------------