

PFA Cable Tidy IEC Exposed Junction Thermocouple with Fitted Mini Plug - Type K -PFA Flat Pair - 2M Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.



Cable Tidy Thermocouple with Fitted Miniature IEC Plug, PFA Twin Twisted or Flat Pair.

A complete precision fine wire, high sensitivity exposed junction thermocouple, with a fitted miniature connector.

Specifications

Specifications

| Product Code | XE-5603-001 |
|---------------------|--|
| General Description | A complete, precision fine wire, high sensitivity thermocouple with a fitted miniature connector and cable tidy. |
| Sensor Type | Exposed Junction |
| Thermocouple Type | K IEC |
| Cable Length | 2M |
| Cable Type | PFA Flat Pair insulated wire |
| Termination Type | Miniature flat pin plug, colour coded in accordance IEC 584 |
| Max. Temperature | Rated 220°C continuous |
| Tolerance | Class 1 |
| Standards Met | IEC 584 |
| Connector Type | Plug |
| Connector Size | Miniature |