



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

PFA Cable Tidy IEC Exposed Junction Thermocouple with Fitted Mini Plug - Type T - PFA Flat Pair - 2M Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Cable Tidy Thermocouple with Fitted Miniature IEC Plug, PFA Twin Twisted or Flat Pair.

A complete precision fine wire, high sensitivity exposed junction thermocouple, with a fitted miniature connector.

Need Your Temperature Sensor calibrated?

Traceable to ISO17025 (UKAS) calibration standards, Labfacility offer a choice of a 3 or 5 point traceable calibration using our in-house calibration facility, you can select any temperature points between -10°C & +200°C.

Please be aware that if you add a calibration option to your order you will be contacted regarding delivery information

Specifications

Specifications

Product Code	XE-5605-001
General Description	A complete, precision fine wire, high sensitivity thermocouple with a fitted miniature connector and cable tidy.
Sensor Type	Exposed Junction
Thermocouple Type	T IEC
Cable Length	2M
Cable Type	PFA Flat Pair Insulated Wire
Termination Type	Miniature flat pin plug, colour coded in accordance IEC 584
Max. Temperature	Rated 220°C continuous
Tolerance	Class 1
Standards Met	IEC 584
Connector Type	Plug
Connector Size	Miniature