

ANSI PFA Cable Tidy Exposed Junction Thermocouple with Fitted Mini Plug - Type T - 2M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Need Your Temperature Sensor calibrated?

Traceable to ISO17025 (UKAS) calibration standards, Labfacility offer a choice of a 3 or 5 point traceable calibration using our in-house calibration facility.

Please be aware that if you add a calibration option to your order you will be contacted regarding delivery information.

Specifications**Specifications**

Product Code	XE-5609-001
General Description	A complete, precision fine wire, high sensitivity thermocouple with a fitted miniature connector
Sensor Type	Cable Tidy Thermocouple with Fitted Miniature ANSI Plug
Thermocouple Type	T
Thermocouple Junction	Exposed
Cable Length	2M
Cable Type	PFA twin twisted
Cable Termination	Miniature Plug
Max. Temperature	220°C
Tolerance	Class 1 wire for high accuracy
Standards Met	ANSI colour coded
Connector Type	Plug
Connector Size	Miniature
Colour	Blue