

Industrial Thermometer Therma 1 (Type K)



[Click here to view product](#)

Product Code	XE-3988-001
General Description	Housed in a robust ABS case that contains 'Biomaster' additive that helps to reduce bacterial growth.
Type	Industrial Digital Thermometer
Thermocouple Type	Type K
Accuracy	±0.4°C ±0.1%
Input / Ranges	0.1°C -99.9 to 299.9°C. 1°C 300 to 1372°C
Resolution of data display	0.1°C & 1°C
Dimensions	25 x 56 x 128mm
Power Supply	Battery 3 x 1.5 volt AAA. Battery life 10000 hours/min 5 years

Sensor Type	Sensor type K miniature thermocouple
Weight	130 grams