



# **LABFACILITY**

TEMPERATURE & PROCESS TECHNOLOGY

## **IEC Bolt Thermocouple, Glassfibre Stainless Steel Over Braided Cable - Type J - M8x1mm Thread Pitch - 13mm Bolt Length**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Specifications****Specifications**

<b>Product Code</b>	XE-5557-001
<b>General Description</b>	Bolt sensors are generally screwed into the surface which normally increases the speed of response and sensitivity of the measurement. Suitable for use up to +250°C
<b>Thermocouple Type</b>	J
<b>Thermocouple Junction</b>	Bolts have grounded junctions
<b>Cable Length</b>	2M
<b>Cable Type</b>	Glassfibre stainless steel over braided cable
<b>Cable Termination</b>	Bare tails
<b>Max. Temperature</b>	250°C
<b>Standards Met</b>	IEC
<b>Bolt thread/pitch</b>	M8 x 1mm