

# GM1311 High Performance Professional Digital Food Thermometer



[Click here to view product](#)

|                     |  |
|---------------------|--|
| <b>Product Code</b> | XE-4002-001  |
| <b>Accuracy</b>     | -10 - 100°C: +/-1°C / 100 - 200°C: +/-1.5°C  |
| <b>Resolution</b>   | 0.1°C / °F   |
| <b>Colour</b>       | White  |
| <b>Dimensions</b>   | 9.29 in x 1.34 in x 0.71 in (23.6 cm x 3.4 cm x 1.8 cm). Display dimension:44x18mm |
| <b>Function #1</b>  | LCD display  |
| <b>Function #2</b>  | Temperature lock   |
| <b>Function #3</b>  | Max/Min measuring  |
| <b>Function #4</b>  | Conversion between Celsius degree and Fahrenheit degree                            |
| <b>Function #5</b>  | Auto power off   |

|                         |                                      |
|-------------------------|--------------------------------------|
| <b>Max. Temperature</b> | 300°C (572°F)                        |
| <b>Min. Temperature</b> | -50°C (-58°F)                        |
| <b>Probe Length</b>     | 149.5 x 6mm                          |
| <b>Weight</b>           | Weight: 0.95 oz (27 g)               |
| <b>Power Supply</b>     | CR2032 button battery x 1 (included) |