

GM1311 High Performance Professional Digital Food Thermometer



[Click here to view product](#)

Product Code	XE-4002-001
Accuracy	-10 - 100°C: +/-1°C / 100 - 200°C: +/-1.5°C
Resolution	0.1°C / °F
Colour	White
Dimensions	9.29 in x 1.34 in x 0.71 in (23.6 cm x 3.4 cm x 1.8 cm). Display dimension:44x18mm
Function #1	LCD display
Function #2	Temperature lock
Function #3	Max/Min measuring
Function #4	Conversion between Celsius degree and Fahrenheit degree
Function #5	Auto power off

Max. Temperature	300°C (572°F)
Min. Temperature	-50°C (-58°F)
Probe Length	149.5 x 6mm
Weight	Weight: 0.95 oz (27 g)
Power Supply	CR2032 button battery x 1 (included)