

GM1311 High Performance Professional Digital Food Thermometer



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

PLEASE NOTE: Lead time will be 7-10 working days when choosing calibration

GM1311 High Performance Professional Digital Food/BBQ/Meat Thermometer

Environment temperature :0 to 40°C (32 to 104°F)

Environment humidity:10-90%RH

Rate Voltage: 3V

1 x Food Thermometer

1 x English user manual

1 x CR2032

Specifications

Specifications

Product Code	XE-4002-001
Accuracy	-10 - 100°C: +/-1°C / 100 - 200°C: +/-1.5°C
Resolution	0.1°C / °F
Colour	White
Dimensions	9.29 in x 1.34 in x 0.71 in (23.6 cm x 3.4 cm x 1.8 cm). Display dimension:44x18mm
Function #1	LCD display
Function #2	Temperature lock
Function #3	Max/Min measuring
Function #4	Conversion between Celsius degree and Fahrenheit degree
Function #5	Auto power off
Max. Temperature	300°C (572°F)
Min. Temperature	-50°C (-58°F)
Probe Length	149.5 x 6mm

Weight	Weight: 0.95 oz (27 g)
Power Supply	CR2032 button battery x 1 (included)