

IEC Velcro Thermocouple, PFA with Miniature Plug - Type T



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

| | |
|----------------------------|--|
| Product Code | XE-5611-001 |
| General Description | Velcro thermocouples can be used as a quick and simple method of measuring the temperature around the surface of a pipe. |
| Thermocouple Type | T |
| Cable Length | 1.5M (0.5M Inside Velcro) |
| Cable Type | PFA |
| Core / Strands | 7/0.2mm |
| Termination Type | Miniature flat pin plug, colour coded in accordance IEC 584 |
| Max. Temperature | 150°C |
| Standards Met | IEC |
| Length | Velcro: 500mm |