



LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

IEC Mineral Insulated Dual Thermocouple with Standard Duplex Type K Plug - 6.0mm x 500mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Need Your Temperature Sensor calibrated?

Traceable to ISO17025 (UKAS) calibration standards, Labfacility offer a choice of a 3 or 5 point traceable calibration using our in-house calibration facility.

Please be aware that if you add a calibration option to your order you will be contacted regarding delivery information.

Specifications**Specifications**

Product Code	XE-2015-001
General Description	310 stainless steel sheath. Sheath can be bent/formed to suit many applications and processes.
Sensor Type	Mineral insulated Type 'K' with Duplex Thermocouple Plug
Thermocouple Type	K
Thermocouple Junction	Insulated hot junction
Probe Diameter	6.0mm
Probe Length	500mm
Termination Type	Standard duplex plug termination (200°C)
Max. Temperature	1100°C
Min. Temperature	-40°C
Standards Met	Conforms to IEC 584 specification
Connector Type	Plug
Connector Size	Standard