



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

General Purpose Pt100 Probe with PFA lead - Type PRT - 4.0mm probe diameter - 150mm probe length - 2M cable length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

General Purpose Pt100 Probe with PFA lead - Type RTD / PRT

Multipurpose Pt100 Probe to IEC751

Class B, 4-wire connection

Rigid 316 stainless steel sheath

Teflon® insulated, screened lead

Sensor temperature rating -50°C to +250°C

Specifications**Specifications**

Product Code	XE-3613-001
General Description	Multipurpose Pt100 Probe to IEC751
Sensor Type	Pt100 (100 Ohms @ 0°C) to IEC 751, Class B, 4 wire
Sheath Material	4.0mm diameter 316 Stainless Steel Sheath
Cable Length	2M
Cable Type	Flexible Teflon® insulated, screened
Termination Type	4-wire tails allows connection to any PT100 instrument 2, 3 or 4 wire
Max. Temperature	250°C
Min. Temperature	-50°C
Standards Met	IEC751
Series	General Purpose Pt100