

## **Stainless Steel Compression Fittings - Metric (M Series) Thread**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

To suit a wide range of mineral insulated and other probe type sensors offered by Labfacility.

Stainless steel compression fittings with Metric (M Series) threads are extensively used for connecting pipes and tubes in various industrial applications.

These fittings feature metric threads, standardized to ensure compatibility with metric piping systems.

Made from stainless steel, they offer superior corrosion resistance, durability, and strength, making them suitable for harsh environments such as chemical processing, pharmaceuticals, and food and beverage industries.

### **Specifications**

## **Specifications**

<b>General Description</b>	Thermocouple installation fittings with Metric Thread
<b>Material</b>	Male in-line stainless steel compression type with nut & olive
<b>Thread</b>	M8, M16 & M20 thread available.
<b>Probe Diameter</b>	Suits 1.0mm, 1.5mm, 2.0mm, 3.0mm, 4.5mm, 6.0mm or 8.0mm probe diameters